JAPANESE

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
Drä upp häl

ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE / INDICATIONS: The AIR STIRRUP PLUS PROVIDES SUPPORT TO PREVENT FURTHER INJURY BY EFFECTIVELY LIMITING ANKLE INVERSION AND EVERSION.

CONTRAINDICATIONS: NA

WARNINGS AND PRECAUTIONS:

If you have any of the following signs, contact your healthcare professional immediately.

- Numbness or tingling in your foot.
• Cold, pale toes. (If you are wearing the Ankle Wrap, loosen or remove it. If the problem persists, contact your healthcare professional.)
• Inability to put weight on your injured foot four days after injury.
• Excessive pain that does not lessen with rest.
• Increase in swelling along with a decrease in mobility.

Step 1: Reduce Pain and Swelling

- Raise your foot.
• Wrap cold pack in a moist cloth and place on your ankle.
• Use 15 minutes at a time.

The Ankle Wrap (included in the Air-Stirrup Plus) will help reduce swelling and should be used for the first few days only. The wrap does not provide support.

- Apply directly on your skin.
• Begin at your toes and wrap up your leg.
• Do not stretch the wrap.
• Repeat several times a day.

Step 2: Apply the Air-Stirrup Ankle Brace, then Walk

- Open brace and lift bottom of aircell until heel straps are fully exposed
• Apply heel straps and adjust heel pad width for a snug fit
• Firmly press heel straps and aircells back in place
• Put on sock or ankle wrap
• Place round edge of heel pad under heel
• Align brace sides with ankle
• Secure straps from bottom to top
• Put on lace shoe or sneaker

Wear the brace for the first 2-3 weeks whenever you are putting weight on your injured ankle. Because it can take months for your ankle to heal completely, wear the Air-Stirrup Ankle Brace during activities.

- Squeeze brace sides together with one hand
• Tighten straps from bottom to top with other hand
• Tighten for comfortable support.

ANKLE BRACE FIT ADJUSTMENT

If you experience pinching or uneven pressure, reposition brace or readjust heel pad and straps. Ankle Brace aircells are preinflated and normally do not require adjustment.

HIGH ALTITUDE: WHEN FLYING, READJUST THE STRAPS TO A COMFORTABLE PRESSURE. AT HIGH ALTITUDES THE AIRCELLS WILL EXPAND BEYOND THEIR OPTIMAL LEVEL AND AIR WILL BE REMOVED FROM THE AIRCELLS.

INSERT INFLATION TUBE ABOUT 1 1/2" INTO VALVE ON TOP OF AIRCELL. REMOVE A SMALL AMOUNT OF AIR BY GENTLY SQUEEZING AIRCELL, OR AD A SMALL AMOUNT OF AIR BY BLOWING INTO TUBE (TOO MUCH AIR WILL REDUCE SUPPORT). SEAL VALVE BY SQUEEZING THE VALVE FLAT, AND, SLOWLY REMOVE TUBE. FOLD AND CLOSE VALVE.

Correctly position your ankle...

Your injured ankle must be correctly positioned with the brace on, before standing or walking. While seated, bend your knee and slowly move your injured foot in toward your body until your foot is flat on the floor.

...then Walk.

Slowly stand. Take small steps, putting the heel of your injured foot on the ground first, then your toes.

Step 3: Exercise

When you begin these exercises, you may have additional soreness - this is normal. Progress through the exercises as you allow. If pain becomes excessive and doesn't subside with rest and ice, contact your healthcare professional immediately (see "Warning Signs" on inside back cover).

Toe Point

Extend your leg and slowly point your toes forward, then slowly point your toes toward ceiling. Alphabet Exercise

Cross your legs with your injured leg on top. Draw the alphabet in capital letters with your leg toe.

Motion Exercises

Each exercise should be started immediately. Each exercise should be done 10 times and repeated 3 times a day, or as recommended by your healthcare professional.

Achilles Stretch

Wrap a towel around the ball of your foot and straighten your leg. Slowly pull the towel toward you. Hold for 30 seconds, then release.

Balance Exercises

After you have successfully completed the Motion Exercises, begin the Balance Exercises. For each exercise you will have to stand, using a chair or wall for balance. Do each exercise 10 times and repeat 3 times a day.

Toe Raise

Stand on both feet, holding a chair or wall for balance. Slowly move up and down on the toes of both feet. Raise your uninjured foot. Slowly move up and down on the toes of your injured foot.

Single Leg Stand

Stand with your injured foot flat on the floor and your other foot raised. Balance for up to 30 seconds.

Care

The entire Air-Stirrup Ankle Brace can be hand washed in lukewarm water with mild soap and air dried.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

LATEX FREE

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THIS DEVICE IS ONLY ONE ELEMENT OF THE OVERALL TREATMENT PROGRAM ADMINISTERED BY A MEDICAL PROFESSIONAL. THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PROFESSIONAL.

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...then Walk. Slowly stand. Take small steps, putting the heel of your injured foot on the ground first, then your toes.

Step 3: Exercise. When you begin these exercises, you may have additional soreness - this is normal. Progress through the exercises as you allow. If pain becomes excessive and doesn't subside with rest and ice, contact your healthcare professional immediately (see "Warning Signs" on inside back cover).

Toe Point. Extend your leg and slowly point your toes forward, then slowly point your toes toward ceiling. Alphabet Exercise. Cross your legs with your injured leg on top. Draw the alphabet in capital letters with your leg toe.

Motion Exercises. Each exercise should be started immediately. Each exercise should be done 10 times and repeated 3 times a day, or as recommended by your healthcare professional.

Achilles Stretch. Wrap a towel around the ball of your foot and straighten your leg. Slowly pull the towel toward you. Hold for 30 seconds, then release.

Balance Exercises. After you have successfully completed the Motion Exercises, begin the Balance Exercises. For each exercise you will have to stand, using a chair or wall for balance. Do each exercise 10 times and repeat 3 times a day.

Toe Raise. Stand on both feet, holding a chair or wall for balance. Slowly move up and down on the toes of both feet. Raise your uninjured foot. Slowly move up and down on the toes of your injured foot.

Single Leg Stand. Stand with your injured foot flat on the floor and your other foot raised. Balance for up to 30 seconds.

Care. The entire Air-Stirrup Ankle Brace can be hand washed in lukewarm water with mild soap and air dried.

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© 2006 DJO, LLC